BEFORE






Blow your nose and wash your hands with soap

How to use the Abbott Panbio™ Rapid Antigen test

SET UP




Place filled test tube in the assembled test tube rack



PREPARATION





0,5"




COLLECTION

Insert swab to reach back of throat


Rub nasal swab against upper part of back of throat. Avoid tongue, teeth and gums

Insert nasal swab 1/2 inch from edge of nostril. Roll 5 times in each nostril and keep the nasal swab in place for 10 seconds (left & right).





Gently pinch the tube to extract the fluid while removing the swab






D Screw the cap onto tube

Break the swab at the break-point where there is an indentation. Leave the swab in the tube


TESTING

* Place the test on a flat surface & do not touch it during the waiting time




Place 5 drops in the rapid response cartridge and wait 15 minutes, no more

Here is an example of a positive test, a negative test and an invalid test







If the result is invalid, redo the nasal swab with a new kit.

* If your rapid test is **NEGATIVE**, you must do a second test 48 hours later ♠ If you are POSITIVE, please call the Info-Covid line: 1-888-662-7482

